

Anti-DYKDDDDK-tag Mouse mAb(HRP-Conjugate)



Catalog Number: 169089

Product name

Ig Type

Anti-DYKDDDDK-tag Mouse mAb(HRP-Conjugate)

IgG2b

Specificity

Applications

Transfected

WB, IP

Antibody description

Dilutions

DYKDDDDK-tag Mouse mAb(HRP-Conjugate)

WB: 1/1000

Preparation

Antigen: Monoclonal antibodies are produced by immunizing mouse with synthetic peptide (DYKDDDDK), coupled to KLH.

Formulation

In PBS(pH 7.4) containing with 1mg/ml BSA 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot analysis of extracts from CHO-K1 cells, untransfected (A) or transfected with different FLAG-fusion proteins (B,C,D), using DYKDDDDK-tag-HRP (HRP-Conjugate) Mouse mAb (1:1000 diluted).