

Anti-Neurofilament Medium Mouse antibody



Catalog Number: 169111

Product name

Anti-Neurofilament Medium Mouse antibody

Specificity

Human

Antibody description

Neurofilament Medium Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human NEFM (AA: 381-443) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide and 0.5% protein stabilizer

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

IgG1


Applications

FC

Validations




Western blot analysis using NEFM mAb against



HEK293 (1) and NEFM (AA


Western blot analysis using NEFM mouse mAb against NTERA-2 (1), SK-N-SH (2), and PC-12 (3)



cell lysate.

Western blot analysis using NEFM mAb against human NEFM recombinant protein. (Expected

Catalog Number: 169111



Flow cytometric analysis of SK-N-SH cells using NEFM mouse mAb (green) and negative control (red).

MW is 31.9 kDa)