

# Anti-Prolactin/PRL Mouse antibody



Catalog Number: 169116

## Product name

Anti-Prolactin/PRL Mouse antibody

## Specificity

Human

## Antibody description

Prolactin/PRL Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant fragment of human PRL (AA: 65-173) expressed in E. Coli.

## Formulation

Purified antibody in PBS with 0.05% sodium azide.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

IgG2b


## Applications

ICC/IF, FC


## Validations



Western blot analysis using PRL mAb against HEK293 (1) and PRL (AA



Immunofluorescence analysis of MCF-7 cells using



PRL mouse mAb (green). Blue


Western blot analysis using PRL mAb against

# Anti-Prolactin/PRL Mouse antibody




Catalog Number: 169116

---



Flow cytometric analysis of Hela cells using PRL mouse mAb (green) and negative control (red).



Immunofluorescence analysis of Hela cells using PRL mouse mAb (green). Blue