

Anti-Dppa2 antibody



Catalog Number: 175003

Product name

Anti-Dppa2 antibody

Ig Type

IgG

Specificity

Mouse

Applications

WB

Antibody description

Rabbit polyclonal antibody to Dppa2

Dilutions

WB: 1:500

Preparation

This antigen of this antibody was peptide

Validations

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C
or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal




Fig1.; Western blot analysis on F9 using anti-DPPA2 polyclonal antibody