

Anti-CSNK1G2 antibody



Catalog Number: 175349

Product name

Anti-CSNK1G2 antibody

Specificity

Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Rabbit

Antibody description

Rabbit polyclonal antibody to CSNK1G2

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human csnk1g2 211-320/415

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; Uterus (Mouse) Lysate at 40 ug; Primary: Anti- CSNK1G2 at 1/1000 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 47 kD; Observed band size: 49 kD