

Anti-FAM13C antibody



Catalog Number: 175415

Product name

Anti-FAM13C antibody

Specificity

Human, Mouse, Rat, Pig, Horse, Rabbit, Sheep

Antibody description

Rabbit polyclonal antibody to FAM13C

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human fam13c1 51-150/585

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4°C.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; HeLa Cell (Human) Lysate at 30 ug; Primary: Anti- FAM13C1 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 66 kD; Observed band size: 63 kD




Fig2: Sample: MCF-7 Cell (Human) Lysate at 40 ug; Primary: Anti-FAM13C1 at 1/300 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 66 kD; Observed band size: 66 kD