

Anti-ND1 antibody

Catalog Number: 175581

Product name

Anti-ND1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to ND1

Preparation

This antigen of this antibody was klh conjugated synthetic peptide derived from human mt-nd1 35-135/318

Formulation

Liquid, 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.

Storage

Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB:1:500-2000

Validations




Fig1: Sample:; Cerebrum (Rat) Lysate at 40 ug; Primary: Anti-MT-ND1 at 1/500 dilution; Secondary: IRDye800CW Goat Anti-Rabbit IgG at 1/20000 dilution; Predicted band size: 36 kD; Observed band size: 35 kD