

# Anti-EPG5 antibody



Catalog Number: 175809

## Product name

Anti-EPG5 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to EPG5

## Preparation

This antigen of this antibody was recombinant protein of epg5

## Formulation

Liquid, PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:1000-1:2000

IHC-P: 1:200-1:500

## Validations




Fig1.; Western blot analysis of Rat Liver Tissue using EPG5 Polyclonal Antibody.. Secondary antibody (catalog# HA1001) was diluted at




Fig2.; Immunohistochemical analysis of paraffin-embedded Mouse Liver Tissue using EPG5 Polyclonal Antibody.