

Anti-AMPD2 antibody



Catalog Number: 175850

Product name

Anti-AMPD2 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to AMPD2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human ampd2. aa range:131-180

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using AMPD2 Polyclonal Antibody