

# Anti-IGLL1 antibody



Catalog Number: 175896

## Product name

Polyclonal

Anti-IGLL1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human

## Applications

## Antibody description

WB

Rabbit polyclonal antibody to IGLL1

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cd179b. aa range:26-75

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality



Fig1.; Western Blot analysis of various cells using CD179b Polyclonal Antibody diluted at 1: 1000