

Anti-KLRC4-KLRK1 antibody



Catalog Number: 175898

Product name

Anti-KLRC4-KLRK1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to KLRC4-KLRK1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human klrk1. aa range:111-160

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of 3T3 cells using CD314 Polyclonal Antibody