

Anti-EPHA2 antibody



Catalog Number: 175964

Product name

Anti-EPHA2 antibody

Specificity

Human, Rat

Antibody description

Rabbit polyclonal antibody to EPHA2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human epha2/3/4. aa range:556-605

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations (kD)




Fig1.; Western Blot analysis of various cells using EphA2/3/4 Polyclonal Antibody