

Anti-TECR antibody



Catalog Number: 176015

Product name

Anti-TECR antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TECR

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gpsn2. aa range:259-308

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

IHC-P: 1:100-1:300

Validations




Fig1.; Western blot analysis of 293T lysis using GPSN2 antibody. Antibody was diluted at 1:2000