

Anti-NCAN antibody



Catalog Number: 176124

Product name

Polyclonal

Anti-NCAN antibody

Ig Type**Specificity**

Rabbit IgG

Human

Applications**Antibody description**

WB

Rabbit polyclonal antibody to NCAN

Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was synthesized peptide derived from neurocan . at aa range: 370-450

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Clonality

Fig1.; Western Blot analysis of Jk cells using Neurocan Polyclonal Antibody