

Anti-OR2T35 antibody



Catalog Number: 176179

Product name

Anti-OR2T35 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR2T35

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or2t2/2t35. aa range:61-110

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations



Fig1::; Western Blot analysis of various cells using Olfactory receptor 2T2/35 Polyclonal Antibody diluted at 1: 500