

Anti-OR4C12 antibody

Catalog Number: 176183

Product name

Anti-OR4C12 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR4C12

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or4c12. aa range:260-309

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations



Fig1.: Western Blot analysis of various cells using Olfactory receptor 4C12 Polyclonal Antibody