

Anti-OR5AK3P antibody



Catalog Number: 176210

Product name

Anti-OR5AK3P antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to OR5AK3P

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or5ak3p. aa range:241-290

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations



Fig1:: Western Blot analysis of various cells using Olfactory receptor 5AK3 Polyclonal Antibody