

Anti-OR5B12 antibody

Catalog Number: 176213

Product name

Anti-OR5B12 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR5B12

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or5b12. aa range:195-244

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1:: Western Blot analysis of A549 cells using Olfactory receptor 5B12 Polyclonal Antibody