

Anti-OR7A10 antibody

Catalog Number: 176227

Product name

Anti-OR7A10 antibody

Specificity

Human, Monkey

Antibody description

Rabbit polyclonal antibody to OR7A10

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or7a10. aa range:260-309

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations

Fig1.; Western Blot analysis of various cells using Olfactory receptor 7A10 Polyclonal Antibody