

Anti-OR8D1 antibody

Catalog Number: 176230

Product name

Anti-OR8D1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OR8D1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human or8d1. aa range:234-283

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations



Fig1.: Western Blot analysis of various cells using Olfactory receptor 8D1 Polyclonal Antibody