

Anti-PEG3 antibody

Catalog Number: 176248

Product name

Anti-PEG3 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PEG3

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human peg3. aa range:1031-1080

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of A549 cells using PEG3 Polyclonal Antibody