

Anti-RAP2C antibody

Catalog Number: 176283

Product name

Polyclonal

Anti-RAP2C antibody

Ig Type**Specificity**

Rabbit IgG

Human, Mouse

Applications**Antibody description**

WB

Rabbit polyclonal antibody to RAP2C

Dilutions**Preparation**

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rap2c. aa range:97-146

Validations**Formulation**

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


Clonality

Fig1.; Western Blot analysis of various cells using Rap 2C Polyclonal Antibody diluted at 1: 500