

Anti-TMBIM1 antibody



Catalog Number: 176285

Product name

Anti-TMBIM1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TMBIM1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tmbim1. aa range:251-300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1.: Western Blot analysis of various cells using RECS1 Polyclonal Antibody diluted at 1: 1000