

Anti-RIMS4 antibody



Catalog Number: 176294

Product name

Polyclonal

Anti-RIMS4 antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to RIMS4

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rims4. aa range:33-82

Validations

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality



Fig1.; Western Blot analysis of various cells using Rim4 Polyclonal Antibody diluted at 1: 500