

Anti-RNF130 antibody

Catalog Number: 176295

Product name

Anti-RNF130 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to RNF130

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rnf130. aa range:231-280

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000


Validations

Fig1.; Western Blot analysis of various cells using RNF130 Polyclonal Antibody diluted at 1: 1000