

Anti-RASL10A antibody



Catalog Number: 176296

Product name

Anti-RASL10A antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to RASL10A

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rasl10a. aa range:107-156

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of MOUSE-SPLEEN cells using RRP22 Polyclonal Antibody diluted at 1: 1000