

Anti-RSAD1 antibody



Catalog Number: 176297

Product name

Anti-RSAD1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to RSAD1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human rsad1. aa range:298-347

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.: Western Blot analysis of various cells using RSAD1 Polyclonal Antibody