

# Anti-S100A7L2 antibody



Catalog Number: 176298

## Product name

Anti-S100A7L2 antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to S100A7L2

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human s100a7l2. aa range:38-87

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations




Fig1.; Western Blot analysis of HuvEc cells using S-100A7L2 Polyclonal Antibody