

Anti-SCN4B antibody

Catalog Number: 176302

Product name

Anti-SCN4B antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to SCN4B

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human scn4b. aa range:61-110

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1.: Western Blot analysis of various cells using Scn4b Polyclonal Antibody