

Anti-SERPINB12 antibody



Catalog Number: 176307

Product name

Anti-SERPINB12 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to SERPINB12

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human serpinb12. aa range:241-290

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using Serpin B12 Polyclonal Antibody