

Anti-SCAF4 antibody



Catalog Number: 176308

Product name

Anti-SCAF4 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to SCAF4

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human sfrs15. aa range:601-650

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of COLO cells using SFRS15 Polyclonal Antibody diluted at 1: 2000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA).