

# Anti-HTR3D antibody



Catalog Number: 176319

## Product name

Anti-HTR3D antibody

## Specificity

Human

## Antibody description

Rabbit polyclonal antibody to HTR3D

## Preparation

This antigen of this antibody was synthesized peptide derived from sr-3d . at aa range: 10-90

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IF

## Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

## Validations




Fig1:: Western Blot analysis of HT-29 cells using SR-3D Polyclonal Antibody