

Anti-HTR3D antibody



Catalog Number: 176319

Product name

Anti-HTR3D antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to HTR3D

Preparation

This antigen of this antibody was synthesized peptide derived from sr-3d . at aa range: 10-90

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of HT-29 cells using SR-3D Polyclonal Antibody