

Anti-TAAR3P antibody



Catalog Number: 176332

Product name

Anti-TAAR3P antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to TAAR3P

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human taar3. aa range:206-255

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IF

Dilutions

WB: 1:500-1:2000

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using TAAR3 Polyclonal Antibody