

Anti-TBPL2 antibody



Catalog Number: 176335

Product name

Polyclonal

Anti-TBPL2 antibody

Ig Type

Specificity

IgG

Human

Applications

Antibody description

WB, IHC-P

Rabbit polyclonal antibody to TBPL2

Dilutions

Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tbpl2. aa range:81-130

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Validations

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality




Fig1:; Western Blot analysis of HeLa cells using TBPL2 Polyclonal Antibody