

Anti-TCEAL1 antibody



Catalog Number: 176337

Product name

Anti-TCEAL1 antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TCEAL1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tceal1. aa range:91-140

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P:1:50-300

Validations




Fig1.; Western Blot analysis of various cells using TCEAL1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).