

Anti-TCEAL6 antibody



Catalog Number: 176338

Product name

Anti-TCEAL6 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to TCEAL6

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tceal6. aa range:61-110

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using TCEAL6 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventibiotech, MN, USA).