

Anti-TNS1 antibody

Catalog Number: 176342

Product name

Anti-TNS1 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to TNS1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tensin-1. aa range:1273-1322

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of K562 cells using Tensin-1 Polyclonal Antibody