

Anti-GTF3C2 antibody



Catalog Number: 176345

Product name

Anti-GTF3C2 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to GTF3C2

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tf3c2. aa range:71-120

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1.; Western Blot analysis of various cells using TFIIC110 Polyclonal Antibody