

Anti-OXER1 antibody



Catalog Number: 176346

Product name

Anti-OXER1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to OXER1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human oxer1. aa range:212-261

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P, IF

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

IF: 1:100-1:500

Validations




Fig1.; Western Blot analysis of various cells using TG1019 Polyclonal Antibody