

# Anti-Th antibody



Catalog Number: 176347

## Product name

Polyclonal

Anti-Th antibody

## Ig Type

## Specificity

IgG

Rat

## Applications

## Antibody description

WB, IHC-P

Rabbit polyclonal antibody to Th

## Dilutions

## Preparation

WB: 1:500-1:2000

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tyrosine hydroxylase. aa range:5-54

IHC-P:1:50-300

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Validations

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality




Fig1:: Western Blot analysis of various cells using TH Polyclonal Antibody