

# Anti-TNFAIP1 antibody



Catalog Number: 176355

## Product name

Anti-TNFAIP1 antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to TNFAIP1

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tnaf1. aa range:71-120

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB, IHC-P

## Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

## Validations




Fig1.; Western Blot analysis of various cells using TNF-IP 1 Polyclonal Antibody diluted at 1: 500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000