

Anti-TNFAIP1 antibody



Catalog Number: 176355

Product name

Anti-TNFAIP1 antibody

Ig Type

IgG

Specificity

Human, Mouse, Rat

Applications

WB, IHC-P

Antibody description

Rabbit polyclonal antibody to TNFAIP1

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tnaf1. aa range:71-120

IHC-P: 1:100-1:300

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1:: Western Blot analysis of various cells using TNF-IP 1 Polyclonal Antibody diluted at 1: 500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000