

Anti-TSN antibody



Catalog Number: 176358

Product name

Anti-TSN antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TSN

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tsn. aa range:101-150

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations



Fig1.; Western Blot analysis of various cells using Translin Polyclonal Antibody diluted at 1: 2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000