

# Anti-CACNA1H antibody



Catalog Number: 176362

## Product name

Anti-CACNA1H antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Rabbit polyclonal antibody to CACNA1H

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human cacna1h. aa range:462-511

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

IgG

## Applications

WB

## Dilutions

WB: 1:500-1:2000

## Validations



Fig1:; Western Blot analysis of various cells using T-type Ca<sup>++</sup> CP  $\alpha$ 1H Polyclonal Antibody. Secondary antibody (catalog#: HA1001) was diluted at 1:20000