

Anti-TUBA3D antibody



Catalog Number: 176363

Product name

Anti-TUBA3D antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to TUBA3D

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tuba3c/e. aa range:201-250

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of various cells using Tubulin α -3C/D/E Polyclonal Antibody diluted at 1:2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000