

# Anti-TUBA3D antibody



Catalog Number: 176363

## Product name

Anti-TUBA3D antibody

## Ig Type

IgG

## Specificity

Human, Mouse, Rat

## Applications

WB

## Antibody description

Rabbit polyclonal antibody to TUBA3D

## Dilutions

WB: 1:500-1:2000

## Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human tuba3c/e. aa range:201-250

## Validations

## Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

## Clonality

Polyclonal




Fig1:: Western Blot analysis of various cells using Tubulin  $\alpha$ -3C/D/E Polyclonal Antibody diluted at 1:2000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000