

Anti-ZNF771 antibody



Catalog Number: 176395

Product name

Anti-ZNF771 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to ZNF771

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human znf771. aa range:161-210

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations



Fig1.; Western Blot analysis of various cells using ZNF771 Polyclonal Antibody. Secondary antibody [catalog#: HA1001) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA).