

Anti-PSG1 antibody



Catalog Number: 176419

Product name

Anti-PSG1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PSG1

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the internal region of human psg1. aa range:221-270

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-300

Validations



Fig1:: Western Blot analysis of K562 cells using PSG1 Polyclonal Antibody.. Secondary antibody [catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100