

Anti-INHBE antibody



Catalog Number: 176424

Product name

Anti-INHBE antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to INHBE

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human inhbe. aa range:301-350

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-1:300

Validations




Fig1:: Western Blot analysis of L929 cells using Inhibin β -E Polyclonal Antibody. Antibody was diluted at 1:500. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100